

Automatic Stainless Steel Distiller

DZ5Z / DZ10Z / DZ20Z

Performance and characteristic:

- All of products are made from stainless steel.
- Condenser is made of seamless stainless-steel pipes and so with very good cooling effects & great capacity of water produce.
- The cooled water returning to the vaporizing drum can save energy and water.
- Drain valve can push the water in container, prevent bilge.
- When lack water, will auto stop heating, when the return water reach the normal level, then continue heating.
- Immerge electric heating tube made from copper, heat efficiency is better.



DZ5Z

DZ10Z

DZ20Z

Parameters

Model	Electric Supply	Power	Water discharge	Product size (Paper)	Gross Weight (kg)	Net Weight (kg)
DZ5Z	220V/50HZ	4.5kW	≥5L/h	35X34X78	7	5
DZ10Z	380V/50HZ	7.5kW	≥10L/h	42X38X90	10	8
DZ20Z	380V/50HZ	15kW	≥20L/h	51X36X80	13.5	11.5