

Common Stainless Steel Distiller

DZ5 / DZ10 / DZ20

Performance and characteristic:

- All of products are made from stainless steel.
- Condenser is made of seamless stainless-steel pipes and so with very good cooling effects & great capacity of water produce.
- The cooled water returning to the vaporizing drum can save energy and water.
- Drain valve can push the water in container, prevent bilge.
- Immerge electric heating tube made from copper, heat efficiency is better.



DZ5

DZ10

DZ20

Parameters

Model	Electric Supply	Power	Water discharge	Product size (Paper)	Gross Weight (kg)	Net Weight (kg)
DZ5	220V/50HZ	4.5kW	≥5L/h	35 X 26 X 78	7	5
DZ10	380V/50HZ	7.5kW	≥10L/h	42 X 30 X 90	10	8
DZ20	380V/50HZ	15kW	≥20L/h	51 X 36 X 80	13.5	11.5