

Benchtop Low Speed Refrigerated Centrifuge TD4-CXC / iTD4-CXC

Features

1. **Microprocessor control**, speed control by frequency-conversion and CFC compressor.
2. **Maintenance free** Frequency-conversion motor without carbon dust pollution.
3. **Safety Interlock system** that prevents opening the centrifuge lid during running.
4. **10/99 programs can be stored**, quickly recall the last used program.
5. **Imbalance sensor**, Over speed and over temperature cut off
6. **Easy operation**, LCD touch screen, compressor control.
7. **With Emergency lid lock release**, in case of power off.



Parameters

Model	TD4-CXC (LCD)		iTD4-CXC (Touch Display)	
Max Speed/RCF	5500r/min	5220xg	Timer Range	0~99min
Speed Accuracy	±10r/min		Temp Range	-20°C~ +40°C
Display	LCD display		Max Capacity	4×250ml
Noise	≤60dB(A)		Power Supply	AC220v&110V 50Hz 18A
Weight	80kg		Dimension	620x690x415mm

Rotor:

No.1 AR		Max Speed:5500r/min Capacity:12×10/15ml Max RCF:5220xg	No.2 SR		Max Speed:4200r/min Capacity:24×10/15ml Max RCF:3120xg
No.3 SR		Max Speed:4200r/min Capacity:32×10/15ml Max RCF:3120xg	No.4 SR		Max Speed:4200r/min Capacity:8×50ml Max RCF:3120xg
No.5 SR		Max Speed:4200r/min Capacity:72×2/5ml Max RCF:3020xg	No.6 SR		Max Speed:4000r/min Capacity:48×2/5ml Max RCF:3680xg
No.7 SR		Max Speed:4000r/min Capacity:4/8×100ml Max RCF:3120xg	No.8 SR		Max Speed:4000r/min Capacity:4×250ml Max RCF:3300xg
No.9 MTP		Max Speed:4000r/min Capacity:2x2x96holes Max RCF:2300xg			Note:Rotors can customize made