

DispensMate Bottle-Top Dispenser

For safer and more efficient way to handle chemicals

Features

1. Excellent chemical resistance, inert and high strength components made of PTFE, FEP, BSG, PP
2. Fully autoclavable at 121°C (only liquid accessible parts)
3. Available in 4 dispensing volume ranges from 0.5mL to 50mL
4. Easy to disassemble and clean for periodic maintenance
5. Supplied with most popular adapters (S40, GL32, GL38, GL25, GL28)

Applications

Bottle-top dispensers offer high precision and accuracy with ISO8655 conformed calibration for use in areas such as pharmacology, biotechnology, healthcare, food and beverage, chemistry, environmental monitoring and academy&research.

- Easy-to-use volume adjustment
- Smooth piston design provides effortless dispensing experience
- Easy-to-use closure cap protects user from exposure to reagent
- Adapters to comfortably fit most laboratory reagent bottles
- Adjustable aspirating tube for different bottle heights



Parameters

Cat. No.	Volume Range	Graduation	Accuracy error		Precision error	
	niL	mL	μL	%	μL	%
7032100101	0.5-5	0.1	±25	±0.5	±10	±0.2
7032100102	1.0-10	0.2	±50	±0.5	±20	±0.2
7032100103	2.5-25	0.5	±125	±0.5	±50	±0.2
7032100104	5.0-50	1.0	±250	±0.5	±100	±0.2

Cat. No.	Descriptions
17400131	GL32, Adapter, 45/32mm
17400132	GL38, Adapter, 45/38mm
17400133	GL25, Adapter, 32/25mm
17400134	GL28, Adapter, 32/28mm
17400135	S40, Adapter, 45/40mm
17400137	Reagent Bottle (Brown, 1L)
17900158	Filling tube (m)
17000802	Discharge tube (m)

