

DispensMate-Pro Premium Bottle-Top Dispenser

Designed Safe, Sturdy & Stable dispensing directly from the bottle in the range of volumes up to 100 mL

Features

1. High chemical compatibility with use of inert materials and PTFE piston
 2. Easy to disassemble and clean for periodic maintenance
 3. Fully autoclavable at 121°C to ensure sterility
 4. Available in 5 dispensing volume ranges from 0.5mL to 100mL
 5. Supplied with most popular adapters (S40, GL32, GL38, GL25, GL28)
- Fast, precise, and reproducible volume locking mechanism on a clear scale
 - Reagent recovery function reduces reagent waste and prevents dripping when not in use
 - 6 adapters to comfortably fit most laboratory reagent bottles
 - Adjustable aspirating tube for different bottle heights



Parameters

Cat. No.	Volume Range	Graduation	Accuracy error		Precision error	
	mL	mL	μL	%	μL	%
7032212001	0.5-5	0.1	±25	±0.5	±10	±0.2
7032212002	1.0-10	0.2	±50	±0.5	±20	±0.2
7032212003	2.5-25	0.5	±125	±0.5	±50	±0.2
7032212004	5.0-50	1.0	±250	±0.5	±100	±0.2
7032212005	10-100	2.0	±500	±0.5	±200	±0.2

Cat. No.	Descriptions
17400131	GL32, Adapter, 45/32mm
17400132	GL38, Adapter, 45/38mm
17400133	GL25, Adapter, 32/25mm
17400134	GL28, Adapter, 32/28mm
17000069	S40, Adapter, 45/40mm
17400037	Reagent Bottle (Brown, 1L)
17400021	Filling tube (m)
17100073	Discharge tube (m)
17100553	GL45,adapter, 45mm