

dTrite Electronic Titrator

The dTrite Digital Burette provides accurate, precise and convenient bottle-top titration, as well as optimum operator safety. It combines both magnetic stirrer and titrator functions, making it easier for the operator to achieve satisfactory results in the field of analytical chemistry, food industry, water analysis etc.

FEATURES

1. High quality motor with excellent precision and accuracy
2. Electronic control provides a fatigue-free operation
3. Remote control panel prevents manual disturbance during operation
4. Built-in magnetic stirrer for a complete and easy operation
5. Large LCD display parameters



Parameters

Volume Range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10μL
Volume Accuracy	R=0.2% CV=0.07%
Velocity	16 stages
Operating Temperature Range	10°C-30°C
Quality Control	DIN EN ISO 8655
Control Type	External control, including the stirrer and the burette program