

DW300、DW600 Basic Power Supply

Introduction

DW series Basic power supply comprises 2 models, including: DW300 and DW600, which can meet the conventional horizontal electrophoresis, small and medium-sized vertical electrophoresis, membrane acetate electrophoresis experiment requirement, the lightweight and compact, simple and practical characteristics, make it popular in for basic teaching.

Features

1. Molding shell, touch keys, microprocessor intelligent control.
2. Automatic memory.
3. Standard timing operation.
4. Automatic detection of no-load, over-load.



Parameters

Model	DW300	DW600
Type of Output	Constant voltage or constant current	Constant voltage or constant current
Output range	5~300V, 1~300mA, maximum 90W	5~600V, 1~200mA, maximum 120W
Increment	IV, 1mA	IV, 1mA
Timing Range	1min~9hr, 59min	1min~9hr, 59min
Display	LED	LED
Output jacks	Two sets of output jacks	Two sets of output jacks
Dimension (LxWxH)	280x237x118(mm)	280x237x118(mm)
Weight	3.0 (kg)	3.0 (kg)